

Bülte_Sequenzprotokoll.ST25
SEQUENCE LISTING

<110> TRANSMIT Gesellschaft für Technologietransfer mbH

<120> Tierartsspezifischer und quantitativer Nachweis von ZNS-Gewebe in Fleisch und Fleischerzeugnissen

<130> An.159/Bü/Sc/Ab

<140> DE 10361489.3

<141> 2003-12-23

<160> 6

<170> PatentIn version 3.1

<210> 1

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> Vorwärtsprimer RTGcowM56F2a

<222> (1)..(19)

<223> primer for determination BSE-riskmaterial from cow, sheep and goat

<400> 1
acctgcgacc tggagtctt

19

<210> 2

<211> 15

<212> DNA

<213> Artificial Sequence

Bülte_Sequenzprotokoll.ST25

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> Rückwärtsprimer RTGcowM56R2a

<222> (1)..(15)

<223> determination of BSE-riskmaterial from cow, goat and sheep

<400> 2

ctcgcgcatc tgccg

15

<210> 3

<211> 17

<212> DNA

<213> Artificial Sequence

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> TaqManmgb-Sonde OptiR

<222> (1)..(17)

<223> Determination of BSE riskmaterial from cow, sheep and goat

<400> 3

actcgttcgt gccgcgc

17

<210> 4

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> Vorwärtsprimer RTGpigM56F2

<222> (1)..(20)

Bülte_Sequenzprotokoll.ST25

<223> Determination of BSE riskmaterial from pig

<400> 4

gacctgacgtac gtggagtccc

20

<210> 5

<211> 18

<212> DNA

<213> Artificial sequence

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> Rückwärtsprimer RTGpigM56R2

<222> (1)..(18)

<223> Determination of BSE riskmaterial from pig

<400> 5

tggcgctcct cctgctcc

18

<210> 6

<211> 17

<212> DNA

<213> Artificial sequence

<220>

<223> GFAP-gene Exon 5 and Exon 6

<220>

<221> OptiR TaqManmgb-Sonde

<222> (1)..(17)

<223> Determination of BSE riskmaterial from pig

<400> 6

actcgttcgt gccgcgc

17